Glucosefine[®] Pro

Blutzucker-Mess-System

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Operation Manual

Dear Customer,

Thank you for purchasing the Glucofine® Pro Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read the following contents thoroughly and carefully.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the Glucofine® Pro Blood Glucose Monitoring System to easily monitor your blood glucose levels by yourself anywhere, any time.

If you have other questions regarding this product, please contact your place of purchase or your local representative (see last page for contact).

Important Safety Precautions

Please read before use

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. Do **NOT** under any circumstances use the device on newborns or infants.
5. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your doctor to have the results interpreted.
6. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
7. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
8. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
9. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.
10. Proper maintenance and periodically control solution test are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact your place of purchase or your local representative (see last page for contact).

Keep these instructions in a safe place.

Table of Contents

Before You Begin	
Important Information	6
Intended Use and Principle of Action	7
Kit Contents	8
Meter Overview	9
Overview Display / Test Strip	10
Setting the Meter	11 - 13
The Four Measuring Modes	14
Performing a Control Solution Test	15 - 17
Performing a Blood Glucose Test	
Preparation	18
Collecting a Blood Sample	19
Testing Procedure	20 - 21
Alternative Site Testing (AST)	22
Reviewing Test Results in Memory	23 - 24
Downloading Test Results to a Computer or Smartphone	25
Important Information on Maintenance and Storage	
Battery	26 - 27
Care and Maintenance of the Meter	28
Important Information on Test Strips	29
Cleaning the Lancing Device	30
Important Information on Control Solution	30
Expected Test Results	31
Troubleshooting	
Error Messages	32
Potential Problems and Solutions	33 - 35
Specifications	36
Symbol Reference	37

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Please consult the healthcare professional before use.
- The measurement unit used for indicating the concentration of blood or plasma glucose can either have a weight dimension (mg/dL) or a molarity (mmol/L).

The approximate calculation rule for conversion of mg/dL in mmol/L is:

mg/dL	divided by 18	= mmol/L
mmol/L	times 18	= mg/dL

For example: 1) $120 \text{ mg/dL} \div 18 = 6,6 \text{ mmol/L}$
2) $7,2 \text{ mmol/L} \times 18 = 129 \text{ mg/dL}$

Intended Use and Principle of Action

Intended Uses

This system is intended for use outside the body (*in vitro* diagnostic use) by people with diabetes at home and by health care professionals in clinical settings as an aid to monitoring the effectiveness of diabetes control.

It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples (from the finger, palm, forearm and upper arm), and from venous whole blood.

It should not be used for the diagnosis of or screening for diabetes or testing on newborns.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing.

Principle of Action

The system measures the concentration of glucose in blood with an amperometric glucose biosensor embedded in the meter, that measures the volume of electrical signal generated by the reaction of glucose in blood sample with the metabolizing enzyme system (Glucose dehydrogenase (*E. coli*)) implanted on test strips. The signal is translated into glucose concentration and thereafter is displayed on the LCD screen of the meter.


Kit Contents

1. Blood Glucose Meter
2. Operation Manual
3. Functional Carrying Case
4. Warranty Card
5. Blood Glucose Test Strips
6. Logbook
7. 1 x 1.5 V AAA Battery
8. Lancing Device
9. Sterile Lancets for Blood Sampling

Important Note

Please check your Glucofine® Pro kit upon receiving to ensure that the packaging is not damaged and it contains all of the parts shown above. If this is not the case please return the system to the place of purchase for replacement.

Meter Overview


- 1 Test Strip Ejector**
Eject the used strip by pushing up this button.
- 2 Display Screen**
- 3 Down Button (▼)**
- 4 M Button (M)**
Enter the meter memory, store settings and silence a reminder alarm.
- 5 Up Button (▲)**
- 6 Test Strip Slot / Strip Port Comm**
Insert test strip here to turn the meter on for testing.
Connect strip port cable for data transmission.
- 7 Battery Compartment**

Overview Display / Test Strip

Overview Display

- 1 Test Result
- 2 Ketone Warning
- 3 Error Warning
- 4 Memory Mode
- 5 QC Mode
- 6 Measuring Mode
QC - control solution test
AC - before meal
PC - after meal
Gen - any time of day
- 7 Date
- 8 Time
- 9 Day Average
- 10 Alarm Symbol
- 11 Measurement Unit
- 12 Smiley / Low/High Symbol

- 13 Battery Symbol
- 14 Test Strip Symbol
- 15 Blood Drop Symbol
- 16 Code

Test Strip

- 1 Absorbent Hole
- 2 Confirmation Window
- 3 Test Strip Handle
- 4 Contact bars

Apply blood sample here
Check here if enough blood was applied
Hold test strip here

ATTENTION: The side with white marks should face up when inserting the test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.

NOTE: The Glucofine® Pro meter should only be used with Glucofine® Pro Test Strips. Using other test strips with this meter can produce inaccurate results!

Setting the Meter

Before using your meter for the first time or if you change the meter battery, you should update and check your individual settings. After inserting the battery (see page 26) please take the following steps and check if your individual settings are saved.

Entering the Setting Mode

Start with the meter off (no test strip inserted). If necessary push and hold button **M** for 3 seconds. Press and hold ▲ and ▼ at the same time to enter the setting mode.

1. Setting the date

The sequence of the date setting is: YEAR → MONTH → DAY.

When entering the setting mode the year number is flashing. Press ▲ or ▼ to select the correct year. Press **M** to confirm.

Afterwards the month number is flashing. Press ▲ or ▼ to select the correct month. Press **M** to confirm.

Afterwards the date of day is flashing. Press ▲ or ▼ to select the correct date. Press **M** to confirm.

2. Setting the time format

Press ▲ or ▼ to select the desired time format (12h or 24h). Press **M** to confirm.

Setting the Meter

3. Setting the time

With the HOUR flashing, press ▲ or ▼ to select the correct number.
Press **M** to confirm.

With the MINUTE flashing, press ▲ or ▼ to select the correct number.
Press **M** to confirm.

4. Setting the measuring unit

With the measuring unit flashing, press ▲ or ▼ to switch between mg/dL or mmol/L. Press **M** to confirm.

Please ask your doctor or healthcare professional which measuring unit to choose.

5. Setting the low and high target ranges for measuring modes

In the next step you can set the low and high target ranges for the measuring modes "Gen" (any time of day), "AC" (before meal) and "PC" (after meal).

For every measuring mode first "LOW" value then "HIGH" value is flashing. While flashing press ▲ or ▼ until the desired target value appears. Press **M** to confirm.

6. Setting the beep sound On/Off

With "🔊" showing in the left upper display toggle between flashing "On" or "Off" by pressing ▲ or ▼ to set the beep sound on or off. Press **M** to confirm.

Setting the Meter

7. Deleting the memory

With "dEL" and "M" on the display, press **M** to continue with "no" to keep all results in memory.

To delete all the results, press ▲ to select "yes", and then press **M** to confirm deleting all memory records.

8. Setting the reminder alarm for your testing

Your meter has four reminder alarms.

The meter will display "1" with a clock icon and "On" or "OFF" to show its current status. You can toggle between "On" or "OFF" by pressing ▲ or ▼, then press **M** to set alarm time or skip to next reminder alarm. To set an alarm time, press ▲ or ▼ with the hour/minute flashing in sequence to select the correct hour/minute and confirm by pressing **M**.

NOTICE: When the reminder alarm for your testing beeps, press ▲ to switch it off. Otherwise it will beep for 2 minutes and then switch off.

Congratulations! You have completed all settings! The meter will switch off automatically.

NOTE

- You can leave the setting mode anytime by pressing ▲ and ▼ simultaneously.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.
- The settings done so far will be saved in both cases.

The Four Measuring Modes

The meter provides you with four modes for measuring:

1. "Gen": General testing mode (default) at any time of day / without relation to a meal.
Preset target values: 70 - 130 mg/dL (3,9 - 7,2 mmol/L)
2. "AC" - before meal 🍏: Testing on an empty stomach / no meal for the last 8 hours
Preset target values: 70 - 130 mg/dL (3,9 - 7,2 mmol/L)
3. "PC" - after meal 🍷: Testing within 2 hours since last meal
Preset target values: 70 - 180 mg/dL (3,9 - 7,2 mmol/L)
4. "QC" - Quality Check: Testing with control solution

The appropriate measuring mode should be selected with every blood glucose testing to ensure a corresponding marking and evaluation of each test result.

IMPORTANT

- Please consult your doctor to adapt the preset target ranges to your individual ranges.
- Always activate "QC" mode when performing a control solution test to ensure a correct quality check!
- To learn how to activate different measuring modes please see the test procedures with control solution or with blood samples described on the following pages.

Performing a Control Solution Test

The Glucofine® Pro control solution contains a known amount of Glucose (sugar) that reacts with the enzyme on the test strip. It is used to check, if meter and test strips are working together properly and the results are accurate.

When should the control solution test be performed?


- For mandatory quality check of professionals (POCT)
- At least once a week to routinely check the meter and test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- When your blood glucose test results are inconsistent with how you feel, or when you think your results are inaccurate.
- For instruction or if you want to practice the testing process.
- If you drop the meter or think you may have damaged the meter.

Glucofine® Pro Control Solution is not contained in the meter set. It is available as set (normal / high level) with your local distributor or directly with METRADO GmbH (see last page for contact information).

IMPORTANT!

Please read chapter "Important Control Solution Information" on page 30 and the instructions for use delivered with control solution.

Performing a Control Solution Test


To perform a control solution test you will need the Glucofine® Pro blood glucose meter, a blood glucose test strip and control solution.

Important: It is recommended that the control solution test be done at room temperature 20 °C to 25 °C. Make sure your control solution, meter, and test strips are at this specified temperature range before testing.

1. Insert the test strip to turn on the meter
Insert the test strip into the meter. Wait for the meter to display "☞" and "●".
2. Press the M button to mark this test as a control solution test
With "QC" displayed, the meter will store your test result in memory under "QC".


WARNING: Always activate "QC" mode when performing a control solution test to ensure a correct quality check! In addition this avoids to mix up the blood glucose test results with the control solution test results in memory.

3. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out the first drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap.

Performing a Control Solution Test

Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down from 5 to 0. The test result is displayed after 6 seconds.

NOTE: To avoid contaminating the control solution do not directly apply control solution from the vial onto a strip.


4. Read and compare the result

Compare the test result with the target range printed on the test strip vial for the control solution used (“normal” or “high”). It should fall within this range. If not please read the instructions again and repeat the control solution test. If the control solution test is continuously not in that range, do NOT continue with a blood glucose test. Please contact your specialist retailer or the distributor in your country.

If marked as “QC” test (see step 2.) the result is automatically stored as control solution test result.

NOTE: The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.

Performing a Blood Glucose Test: Preparation

To perform a blood glucose test you will need the Glucofine® Pro blood glucose meter, a blood glucose test strip, the lancing device and a sterile lancet.

Please follow the instructions for use in the lancing device insert for collecting a blood sample and also refer to the instruction leaflet delivered with the test strips.

Preparing the Puncture Site

Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the glucose value obtained. Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

Please follow the suggestions below before obtaining a drop of blood:


- Wash and dry your hands before starting.
- Select the puncture site either at a fingertip or another body parts (please see section "Alternative Site Testing" (AST) on how to select the appropriate sites).
- Rub the puncture site for until it feels warm before penetration.

Performing a Blood Glucose Test: Collecting a Blood Sample

Please follow the instructions for use in the lancing device insert for collecting a blood sample.

- **Fingertip testing:**
Press the lancing device's tip firmly against the lower side of your fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.

- **Blood from sites other than the fingertip (AST):**
Replace the lancing device cap with the clear cap for AST.
Please refer to chapter "Alternative Site Testing" on page 22.

Warning! To reduce the chance of infection:

- Always use a new, sterile lancet. Lancets are for single use only.
- Never share a lancet or the lancing device with another person.
- Please be very careful when handling with blood samples. Improper procedures may cause serious hazards to your health.


Performing a Blood Glucose Test: Testing Procedure

To test your blood glucose please follow the steps below:

1. Insert the test strip to turn on the meter
Insert the test strip into the meter. Wait for the meter to display "☰" and "●".
2. Select the appropriate measuring mode by pressing the button ▲.
3. Obtaining a blood sample
Use the pre-set lancing device to puncture the desired site. The size of the blood drop should be at least as big as ●, which is 0,5 microliter (μL) of volume.

4. Apply the sample
Gently apply the drop of blood to the absorbent hole of the test strip at a tilted angle. The confirmation window should be completely filled if enough blood has been applied. Do NOT remove your finger until you hear a beep sound. The measuring starts (duration: 6 seconds) and is shown by a countdown from 5 to 0.


Performing a Blood Glucose Test: Testing Procedure

NOTE:

- Do not press the punctured site against the test strip and avoid to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be completely filled with blood before the meter begins to count down. NEVER try to add more blood to the test strip after the drop of blood has moved away. Discard the used test strip and retest with a NEW test strip.
- If you have trouble regarding the testing procedure please contact your health care professional, your local distributor or METRADO GmbH (see last page for contact data).

5. Read your test result

The result of your blood glucose test will appear after 6 seconds when the meter has counted down to 0. The symbol 😊 or ☹️ shows if the result is within the target range. The blood glucose result will be stored in the memory automatically.

6. Eject the used test strip

Eject the test strip by pushing the eject button on the side. Use an appropriate bin to dispose of used test strips. The meter will switch off automatically.


Please follow the instructions in the lancing device insert when removing the lancet.

WARNING: The used lancet and test strip may be biohazardous. Please discard them carefully according to your local regulations.

Alternative Site Testing (AST)

You can test on a variety of locations on your body. We recommend the upper arm, forearm and other areas of your hand as alternative puncture sites.

Important: There are limitations with testing on alternative sites. Please consult your doctor or health care professional before you perform AST.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, physical exercise, or any other event, take a blood sample from your finger only.

We strongly recommend that you perform AST **ONLY at the following times:**

- In a pre-meal or fasting state (more than 2 hours since the last meal)
- Two hours or more after taking insulin
- Two hours or more after exercise

Do **NOT** use AST if

- you think your blood glucose level is low,
- you are unaware of hypoglycemia,
- you are testing for hyperglycemia,
- your test results do not match the way you feel,
- your routine glucose results often fluctuate.

Reviewing Test Results in Memory

The meter stores the 1000 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, start with the meter switched off. If necessary please switch off the meter by holding the **M** button for 3 seconds.

Reviewing Single Test Results

1. With the meter switched off press and release **M** or **▲**.
“**M**” will appear on the display and the first reading you see is the last blood glucose result along with date, time and the measuring mode.
2. Press **▲** or **▼** to recall the test results stored in the meter.
Press and hold **M** to switch off the meter.

Reviewing Average Test Results

1. With the meter switched off press and release **▼** to enter memory mode for average results. Your 7-day average result measured in “Gen” measuring mode will appear on the display.
2. Press **▲** to review 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of “Gen” (general), AC (before meal), then PC (after meal). Press and hold **M** to switch off the meter.

Reviewing Test Results in Memory

NOTE:

- Any time you wish to exit the memory, keep pressing M for 3 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are not included in the day average results.
- When you are using the meter for the first time, “---” is displayed in memory mode. This means that there are no results in memory.

System Messages while reviewing Test Results

Message	Meaning		
Lo	< 20 mg/dL (1.1 mmol/L)		
LOW 😞	20 - 70 mg/dL (1.1 and 3.9 mmol/L)		
😊	AC	PC	Gen
	<i>71 - 129 mg/dL (3.9 - 7.2 mmol/L)</i>	<i>71 - 179 mg/dL (3.9 - 9.9 mmol/L)</i>	<i>71 - 129 mg/dL (3.9 - 7.2 mmol/L)</i>
HIGH 😞	AC	PC	Gen
	<i>130 - 239 mg/dL (7.2 - 13.3 mmol/L)</i>	<i>180 - 239 mg/dL (10 - 13.3 mmol/L)</i>	<i>130 - 239 mg/dL (7.2 - 13.3 mmol/L)</i>
KETONE?	≥ 240 mg/dL (13.3 mmol/L)		
Hi	> 650 mg/dL (36.1 mmol/L)		

Note: Values in rows printed in italics only refer to the manufacturer preset target ranges. They are replaced by your individual target ranges (see page 14).

Downloading Test Results to a Computer or Smartphone

You can transfer your test results to your personal computer or smartphone by using a strip port USB cable or bluetooth dongle and the corresponding software. Further information on the software and the USB cable or bluetooth dongle is available with your local distributor or with METRADO GmbH (see last page for contact), where you can also obtain the cable and the dongle.

Connecting to a Personal Computer

1. Installing the software after obtaining the required cable
To download the software please visit the website www.gluceofine-pro.de.
2. Connecting to a personal computer
Connect the strip port cable to a USB port on your computer. With the meter switched off, insert the other end of the strip port cable to the meter data port. "PC" will appear on the meter display, indicating that the meter is in communication mode.
3. Data transmission
To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.


If you want to transfer data with a bluetooth dongle to your smartphone please refer to the instructions for use in the corresponding App (available for free in Apple App Store (iOS) and Google Play Store (Android)).

WARNING: While the meter is in transmission mode, it will be unable to perform a test.

Battery

Your Glucofine® Pro meter comes with one 1.5V AAA size alkaline battery.

Inserting/Replacing the Battery

To replace the battery, please make sure that the meter is turned off.

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old battery and replace it with one 1,5V AAA size alkaline battery.
3. Close the battery cover. If the battery is inserted correctly, you will hear a “beep” afterwards.

Low Battery Signal

The meter will display one of the messages below to alert you when the meter power is getting low:

1. The “” appears in display:
The meter is functional and the result remains accurate, but it is time to change the battery.
2. The “” appears with “E-b”, “Error” and “low”:
The power is not enough to do a test. Please change the battery immediately.

Battery

NOTE:

- Replacing the battery does not affect the test results stored in the memory.
- Batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Care and Maintenance of the Meter

Cleaning The Meter

1. To clean the meter exterior, wipe it with a cloth moistened with tap water or an alcohol swab with 70% alcohol, then dry the device with a soft dry cloth. Do NOT rinse with water.
2. Do NOT use organic solvents to clean the meter.

Meter Storage

- Storage conditions: -20 °C to 60 °C, below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Meter Disposal

The used meter should be treated as contaminated that may carry a risk of infection during measurement. The battery should be removed and the meter should be disposed in accordance with local regulations.

Important Information on Test Strips

1. Store the vials with your Glucofine® Pro test strips between 2 °C and 30 °C and at a relative humidity of 10 - 85 %. Do not freeze.
2. Store your test strips in their original vial only. Do not transfer to another container.
3. Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
4. After removing a test strip from the vial, immediately close the vial cap tightly.
5. Touch the test strip only with clean and dry hands. Use each test strip immediately after removing it from the vial.
6. Do not use test strips beyond the expiration date. This may cause inaccurate results.
7. Do not bend, cut, or alter a test strip in any way.
8. Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Cleaning the Lancing Device

The lancing device should be cleaned after each use.

1. Wipe all lancing device's exterior surface and release button with an alcohol swab with 70% alcohol,
2. Let the lancing device dry thoroughly.

Users should wash hands thoroughly with soap and water after handling the lancing device.

Important Information on Control Solution

- Use only Glucofine® Pro control solutions with your monitoring system.
- Do not use the control solution beyond the expiration date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test is done at room temperature (20 - 25 °C). Make sure your control solution, meter and test strips are at this specified temperature range before testing.
- Shake the vial with control solution well before use, discard the first drop of control solution and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- After inserting the test strip, press the "M" button once to activate the QC mode.
- Store the control solution tightly closed at a temperature between 2 °C and 30 °C. Do NOT freeze.

Expected Test Results

The meter provides you with plasma equivalent results.

Time of day	Normal glucose ranges for people without diabetes
Fasting and before meal	< 100 mg/dL (5,6 mmol/L)
2 hours after meal	< 140 mg/dL (7,8 mmol/L)

Source: American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.

Please consult your doctor to determine a target range that works best for you.

Troubleshooting: Error Messages

Message	What it means	Troubleshooting
E-b	Appears when the battery is too low.	Replace the battery immediately.
E-U	Appears when a used test strip is inserted.	Repeat with a new test strip.
E-t	Appears when ambient temperature is above or below system operation range.	System operation range is 8°C to 45°C (46.4°F to 113°F). Repeat the test after the meter and test strip are in the above temperature range.
E-0 E-A E-E	Problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance
E-C	Appears when the coding errors.	
E-F	Appears when test strip is removed while counting down, or insufficient blood volume.	Review the instructions and repeat test with a new strip. If the problem persists, please contact METRADO GmbH for help.

Potential Problems and Solutions

Problem: The meter does not display a message after inserting a test strip.

Possible cause	What to do
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

Potential Problems and Solutions

Problem: The test does not start after applying the blood sample.

Possible cause	What to do
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing "●" appears on the display.
Defective meter.	Please contact customer services.

Potential Problems and Solutions

Problem: The control solution test result is out of range.

Possible cause	What to do
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature 20°C to 25°C before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.
Improper working of meter and test strip.	

Specifications

Model	Glucosefine® Pro
Dimension & Weight	90.3 (L) x 52.3 (W) x 18 (H) mm, 60 g
Power Source	1 x 1.5V AAA alkaline battery
Display	LCD with backlight
Memory	1,000 test results with respective date and time
External Output	strip port cable communication

- Auto sample loading detection
- Auto test strip insertion detection
- Auto reaction time count-down
- Auto switch-off after 3 minutes without action
- Temperature Warning

Operating Condition	8 °C - 45 °C, relative humidity < 85%
Meter: Storage/ Transportation Conditions	-20 °C - 60 °C, relative humidity < 95%
Test Strip: Storage/ Transportation Conditions	2 °C - 30 °C, relative humidity 10 % - 85%
Measurement Units	mg/dL or mmol/L
Measurement Range	20 - 650 mg/dL (1,1 - 36,1 mmol/L)
Expected service life	5 years

Symbol Reference

Symbol	Explanation	Symbol	Explanation
	Manufacturer		Use by
	CExxxx combined with MD Symbol the product fulfills the requirements of Directive 93/42/EEC on medical devices		Combined with IVD Symbol the product fulfills the requirements of Directive 98/79/EC on in vitro diagnostics medical devices
	Medical Device		In vitro diagnostic medical device
	Consult instructions for use!		Attention! See further instructions for use
	This product meets the requirements of the Regulation (EU) 2017/745		Discard 6 months after opening
	Batch code		Model number
	Serial number		For single use only
	EU authorized representative for Medical devices		For self-testing
	Storage/transport at (humidity limitation)		Storage/transport at (temperature limitation)
	Keep away from sunlight		Discard batteries and electronic equipment according to local regulations. Do not dispose with household waste.
	Importer		Distributor

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6.

Glucose[®]fine Pro

Blood Glucose Monitoring System

REF 4239006 PZN 11537111

Blood Glucose Meter

TD-4183

CE 0123 IVD
Directive 98/79/EC

Test Strips

TD-4363

CE 0123 IVD
Directive 98/79/EC

Lancing Device

TD-5010

CE MD
Regulation (EU) 2017/45

Sterile Lancets

TD-5084

CE 0123 MD
Directive 93/42/EEC

 TaiDoc Technology Corporation
B1-7F., No.127, Wugong 2nd Road, Wugu Dist.,
24888 New Taipei City, Taiwan
www.taidoc.com

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EC REP MedNet EC-REP GmbH