Glucoseofine® Pro

Kontroll-Lösung

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Control solution

Intended use

The Glucoseofine® Pro control solution is only for use with the Glucoseofine® Pro blood glucose monitoring system. It is used to check, if meter and test strips are working together properly. It is important to test with control solution routinely to make sure your results are accurate.

Summary

The Glucoseofine® Pro control solution contains a known amount of sugar, that reacts with the test strip.

A test with control solution is similar to that of a blood test. After inserting the test strip, you only need to press the "M" button on the meter to activate QC mode and then use the control solution instead of a drop of blood. The control solution result should fall within the appropriate range printed on the side of the test strip bottle.

IMPORTANT

The control range may change with each new vial of test strips. Always use the control range on the label of your current vial of test strips.

When should the control solution test be performed?

- For mandatory quality check of professionals (POCT)
- At least once a week to routinely check the meter and test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- When your blood glucose test results are inconsistent with how you feel, or when you think your results are inaccurate.
- For instruction or if you want to practice the testing process.
- If you drop the meter or think you may have damaged the meter.

CAUTION

The control solution range printed on the test strip vial is for control solution only. It is used to test meter and test strip performance. It is not a recommended range for your blood glucose level.

Chemistry

Blood Glucose Control Solution is a red solution that contains less than 0.3 % D-glucose. The solution does not contain any human-derived or biological materials.

Ingredient	Percent
Active ingredient: D-glucose	0.03 % - 0.30 %
Inert ingredients	> 99.7 %

Storage and handling

- For in vitro diagnostic use only.
- Glucoseofine® Pro control solutions was especially developed for the Glucoseofine® Pro blood glucose monitoring system. Please do not use control solution from other brands!
- Check the expiration date on the vial. Do not use if expired.
- Tightly close the control solution and store at room temperature between 2 °C and 30 °C. Do not freeze.
- Use only for 90 days after first opening. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.

Test procedure

To do a test with Glucoseofine® Pro control solution, follow the instructions printed in the operation manual that came with the system.

Accurate results depend upon following instructions carefully. Follow below technical pointers whenever you test with Glucoseofine® Pro control solution:

1. It is recommended that the meter, control solution, and test strips come to room temperature (20 - 25 °C) before testing.
2. Shake the control solution well before using.
3. Discard the first drop of control solution you squeeze out.
4. Wipe off the dispenser tip before reclose to keep it clean.

Expected results

At room temperature, Blood Glucose Control Solution test results should fall within the expected range printed on the test strip vial label.

If test results fall outside the expected range printed on the test strip vial, please check the troubleshooting guide in operation manual and repeat the test.

Results that fall outside the expected range may indicate:

- An error in the way you are doing the test.
- Expired test strips and control solutions.
- A problem with the meter.
- Control solution that is too warm or too cool.
- Failure to shake the control solution vial vigorously.

CAUTION

If the control solution test is continuously not in range, do NOT continue with a blood glucose test. Please contact your specialist retailer or the distributor in your country.

Availability

Glucoseofine® Pro control solution (box with each 4 ml normal / high) is available with your local distributor or directly with METRADO GmbH, Auf Pfuhlst 1, D-66589 Merchweiler.

Our commitment to you

Our goal is to provide you with quality healthcare products and dedicated customer service. If you are not fully satisfied with this product or if you have questions about Glucoseofine® Pro, please contact METRADO GmbH or your specialist retailer.

Symbols

Symbol	Bedeutung	Symbol	Bedeutung
	In vitro diagnostic medical device		Manufacturer
	Consult instructions für use!		Model number
	Storage/transport at (temperature limitation)		Attention! See further instructions for use
	Use by		This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices
	Batch code		Discard 3 months after opening
	For self-testing		Authorized representative in the European Community
	Importer		Distributor

TaiDoc Technology Corporation
B1-7F, No. 127, Wugong 2nd Rd., Wugu Dist.,
24888 New Taipei City, Taiwan
www.taidoc.com

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MedNet EC-REP GmbH