## BLUTZUCKER-TESTSTREIFEN

Version 1.0 2024/09 312-4183100-072

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## BLOOD GLUCOSE TEST STRIPS

### Warnings

- For in vitro diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your operation manual before you use this test strip. Use only Glucosefine® Pro test strips with Glucosefine® Pro meter to obtain accurate results and to be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Collecting blood on alternative sites of the body is not recommended, if the peripheral circulation is impaired, as the showed test result may not match the actual blood glucose level.
- It may apply under the following circumstances: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemia, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep meter set and any equipment away from small children. If swallowed, consult a doctor immediately for advice.

### Intended use

The Glucosefine® Pro test strips, when used together with the Glucosefine® Pro blood glucose monitoring system, allow your blood glucose levels to be measured by yourself at home or by healthcare professionals. It uses freshcapillary whole blood samples from the finger, palm, upper arm and forearm, and from venous whole blood. This system is not intended for use in the diagnosis or screening of diabetes, or testing on newborns.

Professionals may use test strips to test capillary and venous blood sample; home use is limited to capillary whole blood testing.

### Principle of Action

The system measures the concentration of glucose in blood with an amperometric glucose biosensor embedded in the meter, that measures the volume of electrical signal generated by the reaction of glucose in blood sample with the metabolizing enzyme system (Glucose hydrogenase (E. coli)) implanted on test strips. The signal is translated into glucose concentration and thereafter is displayed on the LCD screen of the meter.

### About Alternative Site Testing (AST)

**IMPORTANT:** There are limitations for performing AST. Please read the operation manual and consult the doctor before you perform AST.

Alternative site testing (AST) is when individuals check their blood glucose levels using areas of the body other than the fingertip. We strongly recommend that you perform AST ONLY at the following times:

- During a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do NOT use AST if:

- You think your blood glucose is low.
- You may not notice if you are hypoglycemic.
- Your AST results are inconsistent with the way you feel.
- You are testing for hyperglycemia.
- Your routine glucose results often fluctuate.

To obtain a blood sample from the alternative sites, please rub the puncture site for approximately 20 seconds before following the procedures of "Testing Your Blood Glucose".

### Limitations

- Xylose: Do not test blood glucose during or soon after a xylose absorption test. Xylose in the blood can give falsely elevated results.
- Hematocrit: The hematocrit level is limited to between 20 % to 65 %. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: This test strip must not be used for the testing of newborns.
- Metabolites: Dopamine, L-Dopa, methyl dopa, uric acid, ascorbic acid, and acetaminophen at normal blood concentration do not significantly affect blood glucose test results.
- There are no significant interference in the presence of galactose, maltose, or fructose observed in blood glucose tests.
- Lipemic Effects: Blood triglycerides up to 3000 mg/dL (33.9 mmol/L) do not affect the results significantly, but may affect results at higher levels.
- Altitude Effects: Altitudes up to 10,742 feet (3,275m) do not affect test results.
- Use only heparin for anticoagulation of fresh capillary or venous whole blood.


**The following compounds when determined to be in excess of their limitation may produce elevated glucose results:**

Substance	Limiting concentration (mg/dL)	Therapeutic / Physiologic concentration range (or upper limit) (mg/dL)
Acetaminophen (Paracetamol)	> 6,25	0,45 - 3
Ascorbic Acid (Vitamin C)	> 5	2
Bilirubin (UnConjugated)	> 20	0 - 2
Dopamine	> 1,25	0,03
Levo – Dopa	> 0,7	0,02 - 0,28
Methyl – Dopa	> 1,875	0,1 - 0,5
Reduced Glutathione	> 30	24,25 – 32,2
Tolazamide	2 – 2,5	2 – 2,5
Uric Acid	> 10	2 - 8
Galactose	> 250	< 5
Mannitol	> 5000	0,0128
Mannose	> 125	1,15
Xylose	> 3,125	N/A
Lipemic Samples (Triglycerides)	> 3000	30 - 300
Pralidoxime Iodide	> 5	~ 10 (i.v. Dose 500 mg)


### Storage and handling

**Important: Do not use the test strips if they have expired.**

- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- Store the test strips in a cool, dry place between 2° C and 30° C and below 10 - 85 % relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- Store the test strips in their original vial ONLY. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Use each test strip immediately after taking it out of the vial. Close the vial immediately after taking out a strip.
- Keep the vial closed at all times.
- Do not bend, cut, or alter the test strip.

- Absorbent hole  
Apply a drop of blood here. The blood will automatically be absorbed.
- Confirmation window  
This is where you confirm if enough blood has been drawn into the absorbent hole of the strip.
- Test strip handle  
Hold this part to insert the test strip into the slot.
- Contact bars  
Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

**ATTENTION:**  
The front side of the test strip should face up when inserting the test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.

### Testing your blood glucose

Please wash and dry your hands before performing any testing.

**Step 1**  
Insert the test strip fully into the slot of the meter to power on. Wait until the display shows " " and " ". Select the suitable measuring mode by pressing the button " ".

**Step 2**  
Use the prepared lancing device to puncture the selected site. The sample volume should be about 0,5 µl " ". Touch the blood drop with the absorbent hole of the test strip.  
Do NOT apply a smeared blood sample.  
The blood will immediately be absorbed. As soon as the confirmation window is fully covered, the meter starts testing. Do not remove your finger, until you hear a beep-sound.

**Step 3**  
After 6 seconds (count down from 5 to 0), the meter will display your blood glucose level. The symbol " " or " " shows, if your test result is within your determined target range. The test result will be automatically saved in the meter. Discard the used test strip by actuating the strip ejector into a suitable container to dispose. The meter will automatically turn off after removing test strip.

92  
mg/dL

Please refer to your Owner's Manual for more information.  
The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

### Reading your results

Your blood glucose readings deliver plasma equivalent results and are displayed either in milligrams of glucose per deciliter (mg/dL) or in millimoles of glucose per liter (mmol/L). The measurement range of this meter is 20 - 650 mg/dL (1,1 - 36,1 mmol/L).

### Reference values

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meal	< 100 mg/dL (5,6 mmol/L)
2 hours after meals	< 140 mg/dL (7,8 mmol/L)

Source: American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): p. 16.  
Please consult your doctor to determine a target range that works best for you.

### Questionable or inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
- Make sure the confirmation window of the test strip is completely filled with blood.
  - Check the expiry date of the test strips
  - Check the performance of your meter and test strip with the control solutions.

PLEASE NOTE: Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

### Quality Control Testing with Control Solution

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial or on the package. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

IMPORTANT: The reference range of the control solutions may vary with each new vial or package of test strips. Make sure you check the range on the label of your current vial or on the current package.

### Chemical Components

Glucose dehydrogenase (E. coli)	8 %
Electron shuttle	55 %
Enzyme protector	8 %
Non-reactive ingredients	29 %

### Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only. Professionals may use test strips to test capillary and venous blood sample.

Sample size:	0,5 µL
Reaction time:	6 seconds
System measurement range:	20 mg/dL to 650 mg/dL (1,1 to 36,1 mmol/L)
Hematocrit Range:	20 % to 65 %

### Accuracy

Diabetes experts have suggested that glucose meters should be within ±15 mg/dL (0,83 mmol/L) of the reference method when the glucose concentration is lower than 100 mg/dL (5,55 mmol/L) and be within ±15 % of the reference method when the glucose concentration is 100 mg/dL (5,55 mmol/L) or higher. The tables below display how often the Glucosefine® Pro meter achieves this target. The chart is based on a study carried out on 160 patients (each patient was tested six times which had 960 test results) to see how well the Glucosefine® Pro blood glucose monitoring system performed compared to YSI-2300 reference method results.

**Table 1 Accuracy results for glucose concentration < 100 mg/dL (5,55 mmol/L)**

Within ± 5 mg/dL (Within ± 0,28 mmol/L)	Within ± 10 mg/dL (Within ± 0,55 mmol/L)	Within ± 15 mg/dL* (Within ± 0,83 mmol/L)
68,7 % (206/300)	92 % (276/300)	98,7 % (296/300)

**Table 2 Accuracy results for glucose concentration ≥100 mg/dL (5,55 mmol/L)**

Within ± 5 %	Within ± 10 %	Within ± 15 %*
59,7 % (394/660)	89,1 % (588/660)	98,2 % (648/660)

**Table 3 Accuracy results for glucose concentrations between 30,6 mg/dL (1,7 mmol/L) to 646,5 mg/dL (35,92 mmol/L)**

Within ± 15 mg/dL or ± 15 % (Within ± 0,83 mmol/L or ± 15 %)
98,3 % (944/960)

\*Acceptance criteria in EN ISO 15197:2015, 95 % of all differences in glucose values should be within ±15 mg/dL (0,83 mmol/L) for glucose concentration < 100 mg/dL (5,55 mmol/L), and within ±15 % for glucose concentration ≥ 100 mg/dL (5,55 mmol/L).

Note: When test results are compared to the reference values, difference values below 100 mg/dL (5,55 mmol/L) are expressed in mg/dL or mmol/L, while those above 100 mg/dL (5,55 mmol/L) are in percent.

### User performance

A study evaluating glucose values of blood samples from fingertip and alternative sites (palm, forearm and upper arm) obtained by 158 lay persons showed the following results compared to the YSI 2300:

**Table 1 Difference distribution for glucose concentration < 100 mg/dL**

Tested sites	Difference within ± 5 mg/dL	Difference within ± 10 mg/dL	Difference within ± 15 mg/dL
Fingertip	18/36 (50 %)	34/36 (94,4 %)	36/36 (100 %)
Palm	21/57 (36,8 %)	46/57 (80,7 %)	56/57 (98,2 %)
Forearm	13/57 (22,8 %)	44/57 (77,2 %)	57/57 (100 %)
Upper arm	17/57 (29,8 %)	47/57 (82,5 %)	57/57 (100 %)

**Table 2 Difference distribution for glucose concentration ≥ 100 mg/dL**

Tested sites	Difference within ± 5 %	Difference within ± 10 %	Difference within ± 15 %
Fingertip	56/122 (45,9 %)	101/122 (82,8 %)	122/122 (100 %)
Palm	37/101 (36,6 %)	73/101 (72,3 %)	99/101 (98 %)
Forearm	24/101 (23,8 %)	68/101 (67,3 %)	98/101 (97 %)
Upper arm	20/101 (19,8 %)	62/101 (61,4 %)	96/101 (95 %)

### Precision

The CV (%) is less than 5 % both in intermediate precision and repeatability precision.

### Symbols

Symbol	Explanation
	In vitro diagnostic medical device
	For single use only
	Use by
	Importer
	Authorized representative in the European Community
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices

Symbol	Explanation
	Manufacturer
	Batch code
	For self-testing
	Distributor
	Storage/transport at (humidity limitation)

Symbol	Explanation
	Consult instructions for use!
	Model number
	Discard 6 months after opening
	Storage/transport at (temperature limitation)
	Attention! See further instructions for use

**TaiDoc Technology Corporation**  
B1-7F, No. 127, Wugong 2nd Rd., Wugu Dist.,  
24888 New Taipei City, Taiwan  
www.taidoc.com

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**MedNet EC-REP GmbH**
**METRADO GmbH**